


Vertical Electrophoresis (BL- SCZ9)



Technical specifications:

Gel Size (W×L)	120×145 mm
Gel thickness	1,1.5 (mm)
Comb Wells	16, 23
Gel capacity	2
Buffer requirement	1200 (ml)

Characteristic features:

- Intended for protein purification and analysis, ideal for the dimensional electrophoresis, DGGE and preparative gel electrophoresis.
- Equipped with special processing technology aluminum plates, effectively guaranteed the heat dissipation more even, and holding stable pH during electrophoresis running.
- Non-proliferation electric field influence, guaranteed the electrophoresis strip more clear and straight. Modular double hinged core cassettes, matched with the basing gel caster, ensuring the whole operation process more convenient.
- Retractable connectors compatible with most major power supplies.
- For the temperature control, BL-SCZ8 can be connected with a outer constant temperature circulator.
- Safe operation design, power can be cut off automatically when remove the lid.