

Horizontal Electrophoresis (BL- SPCT)



Technical specifications:

Gel Size (W×L)	60×60mm; 120×60mm; 60×120mm 120×120mm.
Sample volume	3, 6,8,11,13, 18 ,25.
Buffer Volume	550ml
Dimension (L×W×H)	300×160×75 (mm)
Net weight	1.1 (kg)

Characteristic features:

- Abundant buffer not only guarantees cooling effects, but also keeps the pH Value stable during the whole experiment process.
- Unfolded electrodes make repair and parts replacement easy and quick.
- With the sample background color design on the tray, convenient to adding the sample.
- The transparent top lid effectively prevents liquid in the cell from volatilizing and avoids electricity leakage.
- Auto-switch-off when the lid is removed.