

## Horizontal Electrophoresis (BL- SPDT)



### Technical specifications:

Gel Size (W×L)	140×140mm; 140×70mm; 70×70mm.
Sample volume	1, 2, 3, 6, 14, 8,17, 13, 29.
Buffer Volume	900ml
Dimension (L×W×H)	390×187×63 (mm)
Net weight	2.0 (kg)

### Characteristic features:

- Abundant buffer not only guarantees cooling effects, but also keeps the pH Value stable during the whole experiment process.
- Unfolded electrodes make repair and parts replacement easy and quick.
- With the sample background color design on the tray, convenient to adding the sample.
- The transparent top lid effectively prevents liquid in the cell from volatilizing and avoids electricity leakage.
- Auto-switch-off when the lid is removed.