

AUTOCLAVE

Steam Sterilizer Horizontal (Rectangular)

Handsomely Made for Long Service and to provide easier and surer Sterilization for use in Hospitals, Nursing Homes, Pharmaceutical Industries and Laboratories.



TRIPLE SAFETY

The sterilizer is provided with triple safety.

- By a pressure control switch which automatically limits the pressure to the required set valve.
- Two spring loaded safety valves one each at chamber and boiler level.
- A dead weight safety valve at chamber level.
- A PLUG SCREEN prevents the discharge line from choking.
- The plug is easily removable for routing cleaning.

A powerful Ejector (for drying sterilized linen) circulates air through the chamber. The circulating air passes through a corrosion resistant metal wool filter.

A dial thermometer is provided to show the chamber temperature.

AUTOMATIC VACUUM BREAKER is provided to break vacuum case of vacuum formation due to steam condensation.

Operating Temperature and pressure: - Sterilizing Temp.: 121°C.

Sterilizing Pressure: 1.2 to 1.5 kg/Cm², (15 psi to 22 psi)

Specification

Model	Inner dimension	Capacity Ltr.	Heater Load	Sterilization temp & Pressure	Operating Pressure
BL -SSHR - 200	18" x 18" x 36"	200	12 kW	1.2 Kg/cm ² (15-22psi) at 121 °C	From 10 psi to 22 psi
BL -SSHR- 250	18" x 24" x 36"	250	12 kW		
BL -SSHR - 350	24" x 24" x 36"	350	15 kW		
BL -SSHR - 450	24" x 24" x 48"	450	18 kW		
BL -SSHR - 550	24" x 24" x 60"	550	22 kW		
BL -SSHR - 650	24" x 24" x 72"	650	30 kW		
BL -SSHR - 800	30" x 36" x 60"	800	45 kW		

NOTE: CUSTOMIZED SOLUTION ALSO AVAILABLE ON REQUEST

✉ customerservice@biolinkk.com

🌐 www.biolinkk.com

☎ +91111-2555-3658 /4503-5753



Biolinkk

📍 Third Floor, Block B1, Plot
No.25, Community Centre,
Janakpuri, Delhi, Pin 110058 (India)