









SPECIAL OFFER ON PLASTICWARE



PIPPETTE TIPS

	Name	Cat	Description	Size	Price (INR)
	10µL RACKED, FILTERED, STERILIZED	BP-FTR-10-96S	0.1-10µl, clear filter tips, graduated, sterile, rack pack	96 Tips /rack	₹ 3,500.00
	20µL RACKED, FILTERED, STERILIZED	BP-FTR20	1-20µl, Low Binding clear filter tips, graduated, rack, sterile	96 Tips /rack	₹ 3,500.00
	1000µL RACKED, FILTERED, STERILIZED	BP-FTVR-1000S	100-1000ul (extra-long 1250ul design), Low Binding clear filter tips	96 Tips /rack,	₹ 3,500.00



UNIVERSAL FIT PIPET TIPS

	Name	Cat	Description	Size	Price (INR)
	10µL CLEAR TIPS, GRADUATED	BP-TB10-10	0.1-10ul, clear tips, graduated Color: Clear	1000/pkt	₹ 760.00
	200µL CLEAR TIPS, GRADUATED	BP-TB200-10	1-200ul, clear tips, graduated Color: Yellow	1000/pkt	₹ 650.00
	1000µL CLEAR TIPS, GRADUATED	BP-TB1000-10	100-1000ul, clear tips, graduated Color: Blue	1000/pkt	₹ 840.00


PCR TUBES

	Name	Cat	Description	Size	Price (INR)
	0.2 ML PCR TUBES WITH FLAT CAP (FOR PCR APPLICATION)	BL-MB-P02	0.2ml PCR Tubes with Flat Cap	500/pkt	₹ 900.00
	0.2 ML PCR 8-STRIP TUBES WITH ATTACHED CAP	BL-MB-P08A	0.2ml PCR 8-Strip Tubes with Attached Cap	120 strips/ Pkt	₹ 3,890.00



REAL TIME PCR TUBE STRIPS

	Name	Cat	Description	Size	Price (INR)
	0.2ML QPCR 8-STRIP TUBES (WITH OPTICAL CAPS) FOR REAL TIME APPLICATION	BL-MB-Q200	0.2ml qPCR 8-Strip Tubes (With Optical Caps)	120 Strip /pkt	₹ 4,890.00
	0.1ml qPCR 8-Strip Tubes (With Optical Caps)	BL-MB-Q100C	0.1ml qPCR 8-Strip Tubes (With Optical Caps) 0.1ml qPCR 8-Strip Tubes (With Optical Caps) , Low profile For Real Time Application	120 Strips	₹ 2,800.00


REAL TIME SEALING FILM

	Name	Cat	Description	Size	Price (INR)
	MB-QSM QPCR 96 WELL PLATE OPTICAL SEALING MEMBRANE	BL-QSM	qPCR 96 Well Plate Optical Sealing Membrane (pressure-sensitive adhesive film) Compatible with All Real Time Machine	100 P/K	₹ 12,000.00




REAL TIME PCR PLATE

	Name	Cat	Description	Size	Price (INR)
	0.1ML LOW PROFILE QPCR 96 WELL PLATE	BL-MB-Q96-LP	0.1ml Low Profile qPCR 96 Well Plate (sub-semi skirted) Compatible With: <ul style="list-style-type: none"> • ABI 7500 Fast • ABI 7900 Fast • ABI 3730XL • ABI 3500XL • Step One Plus • QuantStudio™ 3/5 Fast • QuantStudio™ Fast 	10 PCS per Pack	₹ 2,000.00
	0.2ML STANDARD PROFILE QPCR 96 WELL PLATE (HEIGHT SKIRTED)	BL-MB-Q96	0.2ml Standard Profile qPCR 96 Well Plate (Height Skirted) Compatible With: <ul style="list-style-type: none"> • ABI 7500 • ABI 7900HT -0.2 ml • QuantStudio™ -0.2ml 	10 PCS per Pack	₹ 2,000.00
	384 WELL QPCR PLATE	BL-B71519	384 Well qPCR Plate Compatible With: <ul style="list-style-type: none"> • ABI 7500 HT • BioRad-CFX • QuantStudio™ 5 • QuantStudio™ 7 • QuantStudio™ 6 Flex • Eppendorf – Master Cycler 	40 Plates per box	₹ 23,500.00
	0.2ml Standard Profile qPCR 96 Well Plate (Semi Skirted)	BL-MB-Q96-S	0.2ml Standard Profile qPCR 96 Well Plate (Semi Skirted) Compatible With: <ul style="list-style-type: none"> • QuantStudio™ 3/5 0.2 ml • QuantStudio™ - 3500XL 	10 Plates /Pack	₹ 2,300.00



PCR PLATE

	Name	Cat	Description	Size	Price (INR)
	0.2 ML 96-WELL PLATE NON-SKIRTED (FOR PCR APPLICATION)	137-174C	0.2ml Plate Non Skirted Natural	50 Plates Pack	₹ 10,000.00

MICROCENTRIFUGE TUBE

	Name	Cat	Description	Size	Price (INR)
	1.5 ML MICROCENTRIFUGE TUBE	BPMCTB015	1.5ml microcentrifuge tubes, clear, graduated	500/unit	₹ 650.00
	0.65 ML MICROCENTRIFUGE TUBE	MCTB005	0.65ml microcentrifuge tubes, clear, graduated	1000/unit	₹ 2,100.00
	2.0 ML MICROCENTRIFUGE TUBE	MCTB020	2.0 ml microcentrifuge tubes, clear, graduated	500/unit	₹ 790.00

CENTRIFUGE TUBES

	Name	Cat	Description	Size	Price (INR)
	15 ML CENTRIFUGE TUBS	LCST15BS	15 ml centrifuge tubes bulk sterile, 2 Case	500 tubes/pkt	₹ 4,323.00
	50 ML CENTRIFUGE TUBES	LSCT50BS	50 ml centrifuge tubes bulk sterile, 2 Case	500 tubes/ Pkt	₹ 6,270.00

**SPECIAL
NEW YEAR
OFFER**

On Plasticware order
worth Rs.23000/- get a
Free Boat Airpod

