

Thermo Shaker (BL-TS-10)

Features

- LCD display: It is easy to setup and use.
- Reliable design with inner extra temperature protection.
- Each thermo shaker is equipped with lifter that is easy to replace different blocks.
- Digital accuracy and convenience.
- Precision wells for uniform thermal transfer.
- Assorted blocks for a wide variety of tube sizes.
- Brushless DC Motor: noiseless long service life maintenance-free.
- Conforms to CE safety standard.
- Accurately control and display time, temperature and speed.
- 12-month warranty.



Specifications

Model	BL-TS-10
Temp control range	RT+5°C~100°C
Time range	1min ~ 99h59min
Temp accuracy	<±0.3°C*
Display accuracy	0.1°C
Temp uniformity	<±0.3°C*
Speed range	200~1800 rpm
Orbit	3 mm
Heating time	<12min (from25°Cto100°C) *
Multi points running	Yes (Max 5 points)
Multi cycle	Yes (Max 99 times)
Auto heating	Yes
Auto running	Yes
Auto resume to run	Yes
Power supply	150W
Voltage	AC220V or AC110V/50-60Hz
Fuse	250V 3A 05x20
Size (mm)	260x195x150 mm
Weight	7.5Kg