


Reciprocal Shaker Water Bath



Features

- Accurate temperature control
- Shaking speed is stable and smooth
- With drainage stopcock for easy cleaning
- PID controlled speed and temperature with LED digital display
- Made with SUS 304 stainless steel body with powder coating
- Tube rack are sold separately for your option

Specification

Model No.	BL-150D	BL-350D
Temp. Range	Amb. +5 ~ 100°C	Amb. +5 ~ 100°C
Temp. Accuracy	± 0.1°C	± 0.1°C
Bath Capacity	22 L	38 L
Shaking Speed	30 ~ 200 rpm	30 ~ 200 rpm
Shaking Stroke	0.5 to 1.5 (mm)	0.5 to 1.5 (mm)
Heater	1200w	1600w
Digital timer	999h 999min 999sec	999h 999min 999sec
Platform Dimension (mm)	300 x 200	400 x 300
Dimension Inner (mm)	500 x 290 x 150	690 x 330 x 310
Dimension Outer (mm)	655 x 350 x 330	830 x 400 x 430
Power	220v 60Hz	220v 60Hz

Optional Accessories

Flask Clamp (ml)	125ml	250ml	500ml	1000ml	2000ml	3000ml	5000ml
Capacity (BL-350D)	24	15	8	6	2	2	1
Capacity (BL-150D)	12	6	4	2	2	1	--
Gabble Cover (supplied)		AC-220 Tube Rack			AC-210 Platform with clamp		
