

ULTRASONIC HOMOGENIZER

BL-JY96-IIN, BL-JY88-IIN, BL-JY92-IIN, BL-JY92-IIDN



Model	BL-JY96-IIN	BL-JY88-IIN	BL-JY92-IIN	BL-JY92-IIDN
Frequency	20-25KHZ	20-25 KHz	20-25KHZ	20-25KHZ
Display	LCD	LCD	LCD	LCD
Adjustable Power Range	150W(5%-99%)	250W(5%-99%)	650W(1 %-99%)	900W(1 %-99%)
Accompany Tip (φ)	6mm	6mm	6mm	6mm
Optional Tip ((φ)	3mm	3mm	2, 3, 10, 12, 15mm	2, 3, 10, 12, 15mm
Processing Volume	10-100ml	10-200ml	0.5-500ml	0.5-600ml
Duty Ratio	1-99%	1-99%	1 -99%	0.1-99,9%
Temperature Alarm	No	No	No	0-99°C (Prevent sample overheating)
Alarm	time, overload	time, overload	time, overload	time, temperature, overload
Power Supply	220/110V 50Hz/60Hz	220/110V 50Hz/60Hz	220/110V 50Hz/60Hz	220/110V 50Hz/60Hz
Size of Power Supply Chassis	430*255*300mm	430*255*300mm	430*255*300mm	430*255*300mm
Weight (The Host and Transducer)	9.3Kg	9.8Kg	14Kg	14.2Kg
Size of Sound Abating Trap	345*340*570mm	345*340*570mm	345*340*570mm	345*340*570mm

ULTRASONIC HOMOGENIZER

BL-JY98-IIIN, BL-JY98-IIDN*, BL-JY99-IIDN



Model	BL-JY98-IIIN	BL-JY98-IIDN	BL-JY99-IIDN
Frequency	19.5-20.5 KHz	19.5-20.5 KHz	19.5-20.5KHz
Display	7-inch TFT touch screen	7 inch TFT touch screen	7 inch TFT touch screen
Adjustable Power Range	1200W(20%-99%)	1200W(20%-99%)	1800W(20%-99%)
Accompany Tip	20mm	20mm	25mm
Optional Tip	15, 25mm	15, 25mm	20, 28mm
Processing Volume	50-1000ml	50-1000ml	50-1200ml
Duty Ratio	1-99%	0.1-99.9%	0.1-99.9%
Temperature Alarm	No	Yes(Prevent sample overheating)	Yes (Prevent sample overheating)
Alarm	Over-temperature, overtime, overload	Over-temperature, overtime, overload	Over-temperature, overtime, overload
Size Of Power Supply Chassis	430*255*300mm	430*255*300mm	430*255*300mm
New Weight	12.1 Kg	12.1 Kg	13Kg
Weight (The Host And Transducer)	14.5Kg	14.5Kg	15Kg
Size Of Sound Abating Trap	345*340*570mm	345*340*570mm	345*340*570mm