

Liquid Nitrogen Container WITH RACKS (YDS-35- 125-FS)

The transport storage series are economical small liquid nitrogen tanks designed specifically for biological samples static storage and for transporting biological samples. The series is made of high strength lightweight aluminum alloy, containing multi-layer insulation to make the product safe, lightweight and efficient, and it has a variety of accessories optional.



Features

- High strength, light weight and small size aluminum structure;
- Equipped with belt;
- Ultra-low evaporation loss;
- Large straw capacity;
- Vials tube is optional;
- Lockable lid is optional to prevent unauthorized open;
- Level monitoring system is optional;
- Roller base is optional;
- Liquid nitrogen pump is optional;
- CE certified;
- Five years vacuum warranty;



YDS-47-27



YDS-50-125



YDS-65-216



YDS-65-216L



YDS-35-125

Specifications

MODEL	YDS-10-125	YDS-35-125	YDS-47-127	YDS-50-125	YDS-65-216	YDS-65-216L
MAXIMUM STORATE CAPACITY						
RACK NUMBER	1	6	6	6	4	6
2.0 ml CRYO VIAL (100 cells boxes)	–	–	–	–	2000	2400
2.0 ml CRYO VIALS (25cells boxes)	100	750	750	900	–	–
Shelf	4	5	5	6	5	4
PERFORMANCE						
Capacity L	10	30	47	50	65	65
Static evaporation rate L/day	0.43	0.37	0.45	0.43	0.78	0.78
Static holding time/d	23	97	105	115	83	83
UNIT DIMENSION						
Neck DIA. (mm)	125	125	127	125	216	216
Height (mm)	535	747	716	771	715	700
Outer DIA (mm)	303	451	508	508	575	678
Empty weight (kg)	6.5	18.7	20	20.3	29	41.6



YDS-175-216



YDS-95-216



YDS-115-216



YDS-120-216



YDS-145-216

MODEL	YDS-95-216	YDS-115-216	YDS-120-216	YDS-145-216	YDS-175-216
MAXIMUM STORATE CAPACITY					
RACK NUMBER	6	6	4	6	6
2.0 ml CRYO VIAL (100 cells boxes)	3000	3600	4000	4800	6000
2.0 ml CRYO VIALS (25cells boxes)	–	–	–	–	–
Shelf	5	6	10	8	10
PERFORMANCE					
Capacity L	95	120	121	145	175
Static evaporation rate L/day	0.81	0.85	0.85	0.85	0.87
Static holding time/d	117	141	142	171	202
UNIT DIMENSION					
Neck DIA. (mm)	216	216	216	216	216
Height (mm)	950	860	1050	1020	1050
Outer DIA (mm)	678	575	678	678	678
Empty weight (kg)	43.8	51.4	38.7	51.8	55.7



+91-11-45035753
+91-11-25553658



customerservice@biolinkk.co
www.biolinkk.com



Biolinkk | # Third Floor, Block B1, Plot No.25,
Community Centre, Janakpuri, Delhi-110058



bio_linkk



biolinkkm



biolinkk



biolinkk



biolinkk.com