

Liquid Nitrogen Container



Introducing the CryoSeal Liquid Nitrogen Container IX-2- a premium-grade 1.75-liter container designed for storing and transferring biological samples. Ideal for laboratories and research facilities, this container ensures top-notch quality and reliability at the best prices in India.

Key Features:

Container Type: Open Top Container
Capacity (litres): 1.75 liters
Material: Metal, Fiber

The CryoSeal Liquid Nitrogen Container IX-2 stands as a benchmark in the realm of cryogenic storage solutions. Crafted meticulously with a 1.75-liter capacity, this container merges cutting-edge design with exceptional functionality. Engineered specifically for the secure transfer of biological samples, it alleviates the challenges of manual handling, ensuring the preservation of delicate specimens.

This open-top container is meticulously constructed with a blend of robust metal and fiber materials, guaranteeing durability and longevity. Its versatile design accommodates various sample types and storage requirements, catering to the diverse needs of laboratories, medical institutions, and research facilities.

Our commitment to quality is unwavering. Each CryoSeal Liquid Nitrogen Container IX-2 undergoes rigorous testing, adhering to stringent industry standards. Rest assured, this container embodies excellence, reliability, and precision, promising optimal performance in safeguarding valuable biological samples.

Experience the pinnacle of cryogenic storage with our CryoSeal Liquid Nitrogen Container IX-2. Elevate your sample preservation capabilities with a superior blend of innovation and quality assurance, exclusively available through our company.

Specifications

Model		IX-2
Capacity	Ltrs	1.75
Static Evaporation Loss Rate	LTRS/PER DAY	0.08
Static Holding Time Working	DAYS	21
Working Time	DAYS	13
Total Height	Mm	390
Overall Diameter	Mm	168
Neck tube diameter	Mm	34.5
Empty Weight	KG's	1.8
Full Weight	KG's	3.2
No. of Canister	PCS	5
Canister OD	Mm	25
Canister Height	Mm	120
Straw Holding Capacity Single level 0.25ml	pcs	150
Straw Holding Capacity Single level 0.5ml	pcs	75
Straw Holding Capacity Double level 0.25ml	pcs	-
Straw Holding Capacity Double level 0.5ml	pcs	-