



ECONOMY PRECISION BALANCE



SAB-E



SAB



Features

- ★ Six Level Anti-Shock Adjustable Function.
- ★ One Touch Key Calibration Function.
- ★ 17 Different Weighing Units : gm, ct, oz, lb, ozt, dwt, & Piece Counting Functions.
- ★ Bright and Large & Clear LCD Display with Back light.
- ★ Digital Filtering for better stability result on every environment.
- ★ Default and Auto Lock Working Weighing Mode.
- ★ Optional : RS-232 Standard, USB Interface.
- ★ Optional : Re-Chargeable Battery Back up.

MODEL	CAPACITY (G)	RESOLUTION (G)	PAN (MM)	REPEATABILITY (G)	LINEARITY (G)	RESPONSE TIME (SEC)
-------	--------------	----------------	----------	-------------------	---------------	---------------------

External Calibration

BL-SAB 10 E	120	0.001	∅ 80	± 0.002	± 0.002	< 3
BL-SAB 20 E	220	0.001	∅ 80	± 0.002	± 0.002	< 3
BL-SAB 30 E	300	0.001	∅ 80	± 0.002	± 0.003	< 3
BL-O " "	610	0.01	∅ 128	± 0.01	± 0.01	< 2
BL-SAB	1100	0.01	∅ 128	± 0.01	± 0.01	< 2
BL-SAB 1	1500	0.01	∅ 128	± 0.02	± 0.03	< 2

* Do not include glass windshield as standard accessories for SAB 602, SAB 1102, & SAB 1502.

Optional Accessories



Auxiliary Display



Windshield



Software



Test Weights

BIOLINKK

Third Floor ,Block B1,Plot no.25,
Community Centre, Janakpuri, New
Delhi- 110058