



Vortex Shaker (BL-VM-20)

www.biolinkk.com


 Innovation to support science....



Vortex Shaker (BL-VM-20)



Vortex Shaker Model BL-VM-20 by Biolinkk is a compact, ruggedly constructed instrument designed to deliver consistent and reliable performance for mixing applications.



Highlights

- ◆ Rugged Construction
- ◆ Stability and Safety
- ◆ Versatile Operation

It offers the convenience of both continuous and touch modes, with variable speed control reaching up to 3000 RPM. Its sturdy metallic base with rubber vacuum feet ensures stability during operation, making it ideal for use in laboratories. The shaker includes versatile attachments to accommodate various sample types, enhancing its functionality and ease of use.



Features

- ◆ Rugged Construction: Durable design with a compact structure for long-lasting performance.
- ◆ Stability and Safety: Heavy metallic base with rubber vacuum feet to prevent movement during operation.
- ◆ Versatile Operation: Offers both continuous and touch modes with adjustable speed control.
- ◆ Include one hand insert, micro tube insert & cup attachment



Specifications

Model	BL-VM-20
Speed Range (RPM)	Up to 3000 RPM
Power Supply	AC 220V, 50 Hz $\pm 10\%$
Attachments Included	Hand insert, microtube insert, and cup attachment



Accessories

(Optional)

Description	Packing
Spare cup attachment	1
Spare one-hand attachment	1
Spare one-hand insert	1
Spare microtube insert	1





Add: Third Floor, Plot No. 25, Community Centre, Block B-1, Janakpuri, Delhi-110058

To order contact us:

Email: customerservice@biolinkk.com

Tel: +91-11-20003-5753 / 25553658

www.biolinkk.com

