



Mini Centrifuge

(BL-MINI-6K)

www.biolinkk.com


 Innovation to support science....



Mini Centrifuge (BL-MINI-6K)

Mini Centrifuge (Fixed speed & Variable speed)



The BL-MINI-6K Mini Centrifuge features a sleek design, versatile tube compatibility, automatic operation, and a convenient built-in timing function.



Highlights

- ◆ Innovative Design
- ◆ Enhanced Centrifugal Force
- ◆ User-Friendly Operation
- ◆ Silent, Durable Performance



The BL-MINI-6K Mini Centrifuge boasts a distinctive and innovative design that's both aesthetically pleasing and expertly crafted. It comes equipped with two rotor types and various tube carriers suitable for 2.0ml, 0.5ml, and 0.2ml tubes, along with 0.2ml x 8 PCR tube strips. Its user-friendly design allows for easy operation—the centrifuge starts or stops automatically with the simple action of opening or closing the cover. Additionally, the Mini Centrifuge series includes a built-in timing function for added convenience.



Features

- ◆ **Innovative Design:** Featuring a fully transparent top cover and multiple rotors, catering to diverse experimental needs.

The Mini Centrifuge features a fully transparent rounded top cover and comes equipped with multiple rotors to cater to diverse experimental requirements.

- ◆ **Enhanced Centrifugal Force:** Larger radius rotor design delivers 1.5 times more force than standard rotors at equivalent speeds.

Its centrifugal force is amplified by 1.5 times through a large radius rotor with eight tube slots, delivering 1.5 times greater force compared to a standard six-hole rotor at equivalent rotation speeds.

- ◆ **User-Friendly Operation:** Automatic 95° cover opening, one-hand operation, and hassle-free rotor replacement.

The upper cover conveniently pops open to 95° automatically, triggered by an outwardly designed start/stop press-button. This design facilitates effortless one-handed operation.

- ◆ **Silent, Durable Performance:** Ultra-low noise, maintenance-free motor ensures long-term reliability and safety.

With a unique rotor buckle design, replacing and handling rotors becomes notably convenient. Additionally, the centrifuge operates with an ultra-low noise, durable motor, ensuring silent and stable performance. The motor is maintenance-free, ensuring a long lifespan and safety during use.





Specifications

| | |
|-------------------|---|
| Model | BL-Mini-6K |
| Speed Range | 6000rpm |
| Centrifugal Force | 3000xg |
| Sample Capacity | 8 x 2.0ml/1.5ml/0.5ml/0.2ml tube (additional 0.5ml and 0.2ml tube carriers) 2 x 8 x 0.2ml PCR tube strip |
| Time Range | 1s - 99min59s |
| Working Noise | ≤ 45 dB |
| Voltage | AC220V/AC110V, 50/60HZ |
| Power | 220V/110V 50-60HZ |
| Dimension (WxDxH) | 176x156x121 (mm) |
| Net Weight | ≤1.5kg |



Add: Third Floor, Plot No. 25, Community Centre, Block B-1, Janakpuri, Delhi-110058

To order contact us:

Email: customerservice@biolinkk.com

Tel: +91-11-4503-5753 / 25553658

www.biolinkk.com

