



Electrophoresis Power Supply (BL-PS115)

www.biolinkk.com


 Innovation to support science....



Electrophoresis Power Supply (BL-PS115)



Highlights

- ◆ Microcomputer-controlled power supply with stable functions and high efficiency.
- ◆ Elegant design with rubber key buttons for easy operation.
- ◆ Provides constant voltage/current with timer control (0-999 minutes).
- ◆ Power failure memory function allows resumption from the point of interruption.

The Biolinkk BL-PS115 Electrophoresis Power Supply is a high-performance, microcomputer processor-controlled device designed for reliable and efficient electrophoresis experiments. It boasts a stable functionality, enhanced security, and high efficiency, with a lightweight and compact design that minimizes heat dissipation. The advanced design prevents overheating under different working conditions, ensuring long-term durability and stable operation. The BL-PS115 is equipped with a user-friendly interface featuring rubber key buttons for comfortable parameter setting. The power supply offers constant voltage or current with timer control ranging from 0 to 999 minutes, automatically stopping the output when the set time is completed. In case of a power failure, it can store the elapsed time and resume the experiment from the interruption point. With its streamlined housing and high-strength materials, the BL-PS115 combines elegance with functionality, making it an ideal choice for everyday electrophoresis experiments.



Features

- ◆ Stable functions with microcomputer processor control.
- ◆ Constant voltage or current operation for electrophoresis.
- ◆ Timer control from 0 to 999 minutes with auto-stop function.
- ◆ Power failure memory for resuming experiments.
- ◆ Streamlined housing design with high-strength materials.
- ◆ Rubber key buttons for comfortable and easy parameter setting.
- ◆ Lightweight and compact with efficient heat dissipation.
- ◆ High performance-price ratio for everyday electrophoresis needs.



Specifications

Model	BL-PS115
The output voltage	10-300V (1V continuously adjustable)
Output voltage accuracy	±1%
Output current	4-400mA (1mA continuously adjustable)
Output current accuracy	±1%
Output Power	75W (maximum)
timing	1-999min (1min adjustable)
show	led
Automatic power-off recovery function	have
Output type	Constant voltage or constant current
size	31x24x8cm (LxWxH)
weight	1.5kg
safety certificate	Comply with EU CE certification standards





Add: Third Floor, Plot No. 25, Community Centre, Block B-1, Janakpuri, Delhi-110058

To order contact us:

Email: customerservice@biolinkk.com

Tel: +91-11-20003-5753 / 25553658

www.biolinkk.com

